



Time	Notes	Freestyle Libre Data (invasive) mg/dL	Know Labs Prototype (non-invasive) Averaged Data	Know Labs Prototype Averaged Data normalized to mg/dL	Know Labs Prototype Raw Data normalized to mg/dL
	11-28 Know Labs CGM Test				
	Ingestion of 32mg Glucose at this time stamp				
11/28/2018 11:23:54		93	-0.912	106	70
11/28/2018 11:30:19		93	-0.968	95	86
11/28/2018 11:36:42		98	-1.043	80	105
11/28/2018 11:43:03		106	-0.918	105	94
11/28/2018 11:49:23		123	-0.799	129	75
11/28/2018 11:55:46		155	-0.7	149	136
11/28/2018 12:02:06		165	-0.626	164	208
11/28/2018 12:08:28		175	-0.572	175	217
11/28/2018 12:14:51		196	-0.423	205	209
11/28/2018 12:21:18		201	-0.423	205	171
11/28/2018 12:27:45		202	-0.497	190	198
11/28/2018 12:34:08		192	-0.574	175	226
11/28/2018 12:40:31		169	-0.62	165	163
11/28/2018 12:46:52		155	-0.663	157	126
11/28/2018 12:53:18		143	-0.759	137	135
11/28/2018 12:59:52		137	-0.852	119	118
11/28/2018 13:06:22		132	-0.813	126	114
11/28/2018 13:12:44		122	-0.873	114	112
11/28/2018 13:19:02		125	-0.897	109	130
11/28/2018 13:25:30		120	-0.889	111	113

Freestyle Libre® vs Know Labs - mg/dL Blood Glucose Measurements




Time	Notes	Freestyle Libre Data (invasive) mg/dL	Know Labs Prototype (non-invasive) Averaged Data	Know Labs Prototype Averaged Data normalized to mg/dL	Know Labs Prototype Raw Data normalized to mg/dL
	11-27 Know Labs CGM Test				
	Ingestion of 32mg Glucose at this time stamp				
11/27/2018 10:09:25		102	-6.522	102	114
11/27/2018 10:13:17		105	-6.445	104	107
11/27/2018 10:17:10		108	-6.087	111	100
11/27/2018 10:21:03		115	-5.985	113	110
11/27/2018 10:24:59		118	-5.722	119	114
11/27/2018 10:28:50		125	-5.292	128	118
11/27/2018 10:32:44		132	-4.985	134	142
11/27/2018 10:36:37		140	-4.801	138	140
11/27/2018 10:40:30		142	-4.597	142	141
11/27/2018 10:44:23		143	-4.645	141	148
11/27/2018 10:48:15		144	-4.743	139	137
11/27/2018 10:52:07		143	-4.691	140	143
11/27/2018 10:55:59		142	-4.82	137	136
11/27/2018 10:59:51		141	-4.79	138	128
11/27/2018 11:03:49		140	-4.869	136	148
11/27/2018 11:07:41		136	-5	134	132
11/27/2018 11:11:34		132	-4.706	140	141
11/27/2018 11:15:26		128	-4.84	137	133
11/27/2018 11:19:21		124	-4.913	135	120
11/27/2018 11:23:18		121	-4.988	134	142
11/27/2018 11:27:11		119	-5.244	129	131
11/27/2018 11:31:02		117	-5.106	132	123
11/27/2018 11:34:56		115	-5.385	126	134
11/27/2018 11:38:50		113	-5.558	122	102
11/27/2018 11:42:43		111	-5.921	115	124
11/27/2018 11:46:35		109	-6.19	109	120
11/27/2018 11:50:28		106	-6.149	110	110
11/27/2018 11:54:23		103	-6.502	103	92
11/27/2018 11:58:16		100	-6.785	97	98
11/27/2018 12:02:08		97	-6.979	93	107
11/27/2018 12:06:01		94	-6.791	97	93
11/27/2018 12:09:53		91	-7.047	91	95

Freestyle Libre® vs Know Labs - mg/dL Blood Glucose Measurements



Time	Notes	Freestyle Libre Data (invasive) mg/dL	Know Labs Prototype (non-invasive) Averaged Data	Know Labs Prototype Averaged Data normalized to mg/dL	Know Labs Prototype Raw Data normalized to mg/dL
11/27/2018 12:13:45		88	-7.298	86	86
11/27/2018 12:17:37		85	-7.6	80	102
Time	Notes	Freestyle Libre Data (invasive) mg/dL	Know Labs Prototype (non-invasive) Averaged Data	Know Labs Prototype Averaged Data normalized to mg/dL	Know Labs Prototype Raw Data normalized to mg/dL
	11-20 Know Labs CGM Test				
	Ingestion of 32mg Glucose at this time stamp				
11/20/2018 15:46:58		93	-0.105	93	90
11/20/2018 15:49:18		96	-0.100	97	107
11/20/2018 15:51:40		99	-0.095	100	105
11/20/2018 15:54:02		104	-0.089	105	88
11/20/2018 15:56:30		107	-0.087	108	102
11/20/2018 15:58:52		110	-0.081	112	92
11/20/2018 16:01:13		113	-0.078	116	124
11/20/2018 16:03:34		116	-0.071	119	135
11/20/2018 16:05:57		119	-0.064	122	127
11/20/2018 16:08:21		125	-0.058	125	119
11/20/2018 16:10:48		131	-0.050	128	133
11/20/2018 16:13:13		137	-0.049	132	140
11/20/2018 16:15:36		139	-0.049	134	131
11/20/2018 16:18:00		141	-0.044	136	125
11/20/2018 16:20:24		143	-0.045	139	125
11/20/2018 16:22:47		144	-0.035	141	151
11/20/2018 16:25:11		146	-0.029	144	154
11/20/2018 16:27:35		147	-0.027	143	152
11/20/2018 16:29:57		148	-0.018	144	147
11/20/2018 16:32:20		147	-0.023	143	140
11/20/2018 16:34:41		145	-0.019	142	123
11/20/2018 16:37:04		143	-0.026	141	132
11/20/2018 16:39:27		141	-0.030	139	143
11/20/2018 16:41:51		139	-0.031	140	150
11/20/2018 16:44:14		137	-0.035	138	144
11/20/2018 16:46:37		134	-0.032	139	133
11/20/2018 16:48:59		131	-0.033	137	128
11/20/2018 16:51:23		128	-0.030	136	138
11/20/2018 16:53:44		126	-0.032	134	143
11/20/2018 16:56:08		124	-0.032	130	151
11/20/2018 16:58:36		122	-0.035	127	131
11/20/2018 17:00:59		119	-0.042	121	120
11/20/2018 17:03:22		116	-0.045	119	100
11/20/2018 17:05:45		112	-0.056	116	108
11/20/2018 17:08:10		109	-0.062	114	110
11/20/2018 17:10:36		108	-0.071	112	114
11/20/2018 17:13:02		107	-0.075	110	107
11/20/2018 17:15:26		106	-0.073	110	110
11/20/2018 17:17:49		107	-0.072	110	114
11/20/2018 17:20:15		108	-0.069	111	116
11/20/2018 17:22:40		110	-0.072	112	110
11/20/2018 17:25:03		111	-0.073	112	104
11/20/2018 17:27:27		112	-0.074	111	103
11/20/2018 17:29:53		113	-0.072	112	114
11/20/2018 17:32:17		114	-0.072	112	125
11/20/2018 17:34:39		115	-0.067	110	128
11/20/2018 17:37:01		116	-0.068	108	113



Time	Notes	Freestyle Libre Data (invasive) mg/dL	Know Labs Prototype (non-invasive) Averaged Data	Know Labs Prototype Averaged Data normalized to mg/dL	Know Labs Prototype Raw Data normalized to mg/dL						
11/20/2018 17:39:25		118	-0.068	108	97						
11/20/2018 17:41:48		120	-0.074	109	105						
11/20/2018 17:44:11		122	-0.079	111	97						
11/20/2018 17:46:35		124	-0.082	115	102						
11/20/2018 17:48:58		126	-0.082	120	118						
11/20/2018 17:51:21		127	-0.073	126	131						
11/20/2018 17:53:44		128	-0.062	132	147						
11/20/2018 17:56:06		130	-0.048	137	144						
11/20/2018 17:58:32		132	-0.035	137	147						
11/20/2018 18:00:54		130	-0.026	137	149						
11/20/2018 18:03:15		128	-0.021	135	154						
11/20/2018 18:05:38		125	-0.021	131	142						
11/20/2018 18:08:03		122	-0.033	125	118						
11/20/2018 18:10:30		119	-0.043	119	103						
11/20/2018 18:12:57		115	-0.055	115	107						
11/20/2018 18:15:24		111	-0.064	111	100						
11/20/2018 18:17:49		108	-0.071	108	104						